

PuroSPIN™ Viral DNA and RNA Purification Kit (96-well format)

- *Viral DNA and RNA Purification from blood, plasma, serum, milk, buccal or nasal swabs, saliva*

#NK056-480

Product Manual

Version 2.4

This product is intended solely for research use only. Full material safety data sheet document is available at www.lunanano.com.

For further information about this and other products please visit our website at www.lunanano.com.

Kit Contents

Component	5 x 96 Preps (#NK056-480)
Lysis Buffer LB-V	200 mL
Wash Buffer WB-R1	300 mL
Wash Buffer WB-R2	300 mL *
Elution Buffer EB	50 mL
Water, Nuclease-free	50 mL
Proteinase K (10 mg/mL), Nuclease-free	6 mL
Carrier RNA	5 x 600 µg
96-well Nucleic Acid Purification Plate	5
96-well Deep Well Lysis Plates (2.2 mL volume per well)	5
96-well Deep Well Wash Plates (1.6 mL volume per well)	5
96-well Elution Plates	5
Sealing Film	25

* After addition of 96-100% Ethanol

Storage

Carrier RNA and **Proteinase K** need to be stored at $\leq -20^{\circ}\text{C}$. Once reconstituted in water, store **Carrier RNA** in one-time-use aliquots at $\leq -20^{\circ}\text{C}$. All other components can be stored at room temperature for at least 1 year.

Additional materials and equipment required

- 96-100% Ethanol
- Pipettes, sterile nuclease-free tips
- 1.5 mL nuclease-free microcentrifuge tubes
- Appropriate vacuum manifold, or centrifuge capable of processing stacked microplates.
- Vortex (if not available, can be replaced with vigorous pipetting)
- Appropriate personal protective equipment (gloves, lab coat)
- Water bath or heating block (heated to 55°C).
- TE Buffer (10 mM Tris-HCl, pH 8.0, 1 mM EDTA) or PBS (137 mM NaCl, 2.7 mM KCl, 10 mM Na_2HPO_4 , 2 mM KH_2PO_4 , pH 7.4)

Kit Description

PuroSPIN™ Viral DNA and RNA Purification Kit is designed for extraction of viral genomic nucleic acids from blood, plasma, serum, buccal or nasal swabs, and saliva. The kit can be used both with RNA and DNA viruses as our buffers formulations and membranes are specifically optimized to bind both RNA and DNA nucleic acids. The kit greatly simplifies the extraction process by incorporating a convenient easy-to-use extraction technology based on nucleic acid binding silica spin-columns. The virus in the sample is first digested and nucleic acid released by the treatment with Proteinase K and **Lysis Buffer LB-V**. Following the lysis step, nucleic acid is loaded onto silica-based membranes and washed to remove protein, salt and other contaminants. Finally, purified nucleic acid is eluted from the membranes into a clean 96 well plate in **Elution Buffer EB**, or **nuclease-free water**. High quality DNA and RNA extracted with the **PuroSPIN™ Viral DNA and RNA Purification Kit** can be used for a range of downstream applications, such as reverse transcription or qPCR. PuroSPIN™ Nucleic Acid Purification Plates have a very high binding capacity and efficiently bind both DNA and RNA nucleic acids.

Important Notes

BEFORE USE:

- Add 96-100% Ethanol to the **Wash Buffer WB-R2**:
 - **#NK056-480** (5 x 96 preps): add **225 mL** of Ethanol
- Reconstitute each tube of 600 µg of lyophilized **Carrier RNA** in 600 µL of **Elution Buffer EB** to make 1 µg/µL solution.
 - Pulse-spin the tube before buffer addition to ensure the pellet is at the bottom.
 - Immediately after reconstitution aliquot the **Carrier RNA** solution into 50 µL aliquots. Store the aliquots at **≤ -20 °C**.
 - When working with Carrier RNA, keep the solution on ice to prevent RNA degradation. Avoid multiple freeze-thaw cycles of the solution aliquots as this might lead to RNA degradation.
 - Reconstitute one tube at a time to prevent **Carrier RNA** degradation



Lysis Buffer LB-V contains **guanidine hydrochloride**. It is a known irritant that is harmful if inhaled, swallowed, or contacted with skin. Wear appropriate protective equipment when working with this reagent.

If precipitate is observed in the **Lysis Buffer LB-V** or **Wash Buffer WB-R1**, re-dissolve it by warming the solution to 37 °C. Cool the solution back to room temperature before use.

RNA is highly prone to degradation by the RNase enzymes. RNases are very stable enzymes that are present on human skin and commonly contaminate lab equipment and surfaces. Therefore, when extracting viral RNA always wear gloves, use RNase-free filtered tips, and pre-treat all non-disposable equipment (pipettes, vortex, centrifuge) and surfaces with RNase removal solution. Ensure all of the buffers are tightly closed when not in use. Take particular care with the **Elution Buffer EB** and **nuclease-free water** solutions, since these are used for storage of eluted RNA.

In addition to possibility of degradation by contaminating RNase enzymes, RNA can also spontaneously break down through auto-hydrolysis. To minimize the probability of both hydrolysis and RNase degradation, the samples should either be kept on ice or frozen prior to RNA extraction. Both bacterial and mammalian cells can be pelleted and stored at -70 °C until the time of RNA extraction. Lysis and homogenization steps should be performed as soon as possible after the sample is removed from ice or thawed. If fresh samples are used, the time between sample collection and RNA extraction should be minimized. Following the extraction, RNA should either be (i) kept on ice and used as soon as possible for the downstream application, or (ii) immediately frozen and stored at -70 °C until needed.

While viral DNA is generally much more stable compared to viral RNA, when extracting the nucleic acid from DNA viruses, it is still recommended to use DNase-free filtered tips and keep sample on ice after extraction.

All sample material and waste should be regarded as potentially infectious. Please wear proper protective equipment when handling samples and waste solutions and follow proper procedures for handling biohazard materials. Work in biosafety cabinet or similar laminar air flow cabinet if possible until samples are lysed. Disinfect all work surfaces thoroughly after the procedure. Follow proper waste disposal guidelines as recommended by local authorities.

Carrier RNA (included with the kit) serves an important role in the purification of viral nucleic acids. First, viral nucleic acids can be present at very low copy number in the sample. At such low concentrations the ability of DNA or RNA to bind to the silica membrane is significantly reduced. However, addition of Carrier RNA enhances the binding of low copy nucleic acids to the membrane, greatly increasing the yield of the purified product. Second, in the case of small contamination of the sample with RNase enzymes, more abundant Carrier RNA preferentially acts as the RNase enzyme substrate, preventing RNase binding and degrading of the viral RNA.

Protocol. Purification of Viral DNA and RNA

Within 30 min of starting the nucleic acid extraction, prepare the following reagents. Keep **Working Buffer LB-V** on ice to prevent **Carrier RNA** degradation.

- Make **Working Buffer LB-V** by:
 - Adding 5 μ L of 1 μ g/ μ L of **Carrier RNA** solution to each 300 μ L of **Lysis Buffer LB-V** used.
 - Adding 10 μ L of 10mg/ml **Proteinase K** solution to each 300 μ L of **Lysis Buffer LB-V** used.
- Calculate the total amount of **Working Buffer LB-V** needed as follows:
 - Volume of **Working Buffer LB-V** needed = (300 μ L) x (number of samples)

A. Extraction Protocol for Centrifuges

Step	Procedure
1	<p>Sample prep:</p> <p>a) <u>Purification from samples in Viral Transport Medium (VTM)</u></p> <ul style="list-style-type: none">▪ Centrifuge VTM tubes containing specimen at maximum speed for 1 min.▪ Use 200 μL collected from the bottom of the tube. <p>b) <u>Purification from blood, serum, plasma</u></p> <ul style="list-style-type: none">▪ Use 200 μL of EDTA or citrate treated blood, serum or plasma. <p>c) <u>Purification from milk</u></p> <ul style="list-style-type: none">▪ Use 200 μL of milk sample. <p>d) <u>Purification from fresh nasal or buccal swabs</u></p> <ul style="list-style-type: none">▪ Swirl the swab for 2-4 minutes in 200 μL of TE or PBS buffer to transfer the cells and viral particles to the buffer. Use the buffer for <u>Step 3</u>. <p>e) <u>Purification from dried nasal or buccal swabs</u></p> <ul style="list-style-type: none">▪ Break the swab tip into the mixture of 300 μL of Working Buffer LB-V and 200 μL of PBS in a well of fresh Deep Well Lysis Plate.▪ Seal the plate with Sealing Film and incubate for 10 min at 55°C.▪ Remove Sealing Film and carefully take out the swab tip and move directly to <u>Step 5</u>. <p>f) <u>Purification from saliva</u></p> <ul style="list-style-type: none">▪ Centrifuge saliva sample at 4,000 x g for 5 min. Discard the supernatant.▪ Resuspend the pelleted cells in 200 μL of TE or PBS buffer. Use the buffer for <u>Step 3</u>.

2	Add 300 μ L of Working Buffer LB-V to each well of one of the Deep Well Lysis Plates .
3	Transfer 200 μ L of sample or resuspension buffer from <u>Step 1</u> to a well of the Deep Well Lysis Plate containing Working Buffer LB-V .
4	Seal the plate with Sealing Film. Mix by vortexing. Incubate the samples for 10 min at 55°C. Remove the Sealing Film.
5	Add 300 μ L of 96-100% Ethanol to each well. Seal with another Sealing Film. Mix by vortexing. Remove the Sealing Film.
6	Transfer 800 μ L of the solution onto the PuroSPIN™ Nucleic Acid Purification Plate placed above a <i>fresh</i> Deep Well Wash Plate . Centrifuge the stacked plate at 5,000 x g for 2 min to load the nucleic acids onto the membrane. Discard the flow-through from the lower plate and reassemble.
7	Add 500 μ L of Wash Buffer WB-R1 to each sample of the plate. Centrifuge at 5,000 x g for 1 min. Discard the flow-through again and reassemble.
8	Centrifuge at 5,000 x g for 1 min to dry out the membrane.
9	Add 500 μ L of Wash Buffer WB-R2 (<i>with Ethanol added</i>) to each sample of the plate. Centrifuge at 5,000 x g for 1 min. Discard the flow-through again and reassemble.
10	Centrifuge at 5,000 x g for 1 min to dry out the membrane.
11	Move the PuroSPIN™ Nucleic Acid Purification Plate onto a fresh 96 well elution plate.
12	Pipette 50 μ L of Elution Buffer EB or nuclease-free water onto the center of each well. Centrifuge at 5,000 x g for 1 min to elute the DNA and RNA. Keep the flow-through that contains the product.

B. Extraction Protocol for Vacuum Manifolds

Individual manifolds vary in their specifications. Check with your manifold manufacturer for specific usage instructions.

Step	Procedure
1	<p>Sample prep:</p> <ul style="list-style-type: none">a) <u>Purification from samples in Viral Transport Medium (VTM)</u><ul style="list-style-type: none">▪ Centrifuge VTM tubes containing specimen at maximum speed for 1 min.▪ Use 200 µL collected from the bottom of the tube.b) <u>Purification from blood, serum, plasma</u><ul style="list-style-type: none">▪ Use 200 µL of EDTA or citrate treated blood, serum or plasma.c) <u>Purification from milk</u><ul style="list-style-type: none">▪ Use 200 µL of milk sample.d) <u>Purification from fresh nasal or buccal swabs</u><ul style="list-style-type: none">▪ Swirl the swab for 2-4 minutes in 200 µL of TE or PBS buffer to transfer the cells and viral particles to the buffer. Use the buffer for <u>Step 3</u>.e) <u>Purification from dried nasal or buccal swabs</u><ul style="list-style-type: none">▪ Break the swab tip into the mixture of 300 µL of Working Buffer LB-V and 200 µL of PBS in a well of fresh Deep Well Lysis Plate.▪ Seal the plate with Sealing Film and incubate for 10 min at 55°C.▪ Remove Sealing Film and carefully take out the swab tip and move directly to <u>Step 5</u>.f) <u>Purification from saliva</u><ul style="list-style-type: none">▪ Centrifuge saliva sample at 4,000 x g for 5 min. Discard the supernatant.▪ Resuspend the pelleted cells in 200 µL of TE or PBS buffer. Use the buffer for <u>Step 3</u>.
2	Add 300 µL of Working Buffer LB-V to each well of one of the Deep Well Lysis Plates .
3	Transfer 200 µL of sample or resuspension buffer from <u>Step 1</u> to a well of the Deep Well Lysis Plate containing Working Buffer LB-V .
4	Seal the plate with Sealing Film. Mix by vortexing. Incubate the samples for 10 min at 55°C. Remove the Sealing Film.

5	Add 300 μ L of 96-100% Ethanol to each well. Seal with another Sealing Film. Mix by vortexing. Remove the Sealing Film.
6a	Assemble your vacuum manifold with the PuroSPIN™ Nucleic Acid Purification Plate placed above the waste tray, or if no waste tray is provided by your manufacturer, a <i>fresh Deep Well Wash Plate</i> can be used in the waste tray's place.
6b	Transfer 800 μ L of the solution onto the PuroSPIN™ Nucleic Acid Purification Plate . Apply vacuum to load the nucleic acids onto the membrane. Aim for 1-2 drops per second. Release vacuum slowly. Discard the flow-through from the lower plate, if needed, and reassemble.
7	Add 500 μ L of Wash Buffer WB-R1 to each sample of the plate. Apply vacuum to wash the membrane. Aim for 1-2 drops per second. Release vacuum slowly. Discard the flow-through from the lower plate, if needed, and reassemble.
8	Apply vacuum for an additional 30 s to dry out the membrane. Release vacuum slowly.
9	Add 500 μ L of Wash Buffer WB-R2 (with Ethanol added) to each sample of the plate. Apply vacuum to wash the membrane. Aim for 1-2 drops per second. Discard the flow-through from the lower plate, if needed, and reassemble. Release vacuum slowly.
10	Apply vacuum for an additional 30 s to dry out the membrane. Release vacuum slowly.
11	Reassemble the vacuum manifold to position the PuroSPIN™ Nucleic Acid Purification Plate above a fresh 96 well elution plate.
12	Pipette 50 μ L of Elution Buffer EB or nuclease-free water onto the center of each well. Briefly apply vacuum to pull the liquid into the membrane. Wait 30 s, then apply vacuum to elute the DNA and RNA. Carefully release the vacuum. Keep the flow-through that contains the product.


Troubleshooting

Problem	Possible causes and solutions
Membrane clogging	<p>Starting material was not completely disrupted. Perform more thorough homogenization. Try a longer digestion with Proteinase K. Use smaller amount of starting material.</p> <p>Precipitates were not completely removed. When using plasma samples, remove visible kryoprecipitates by centrifugation for 3,000 x g for 5 min.</p>
Low nucleic acid yield	<p>Ethanol was not added to the lysate. Make sure the indicated amount of ethanol is added to the sample lysate before loading it onto the membrane.</p> <p>Ethanol was not added to Wash Buffer WB-R2 Add the indicated amount of 96-100%.</p> <p>Too much starting material was used. Reduce the amount of material used for the lysate preparation. Do not use higher amount of material than indicated in the protocol.</p> <p>Carrier RNA was not added to the lysate. Add the indicated amount to create the working buffer before starting each extraction.</p> <p>Degraded Carrier RNA. Avoid multiple cycles of freeze-thawing of reconstituted Carrier RNA. Always keep thawed Carrier RNA on ice. Store frozen aliquots of Carrier RNA below -20 °C</p> <p>Viral nucleic acid eluate too dilute. Use recommended 50 µL of Elution Buffer EB.</p>
Degraded RNA	<p>RNase contamination Use RNase-free tips and tubes. Clean work surfaces and non-disposable items with RNase removal solution. Wear gloves.</p> <p>Inappropriate handling of the sample. When RNA is extracted from fresh samples, the samples should be kept on ice, lysis and homogenization should be performed as soon as possible. Frozen samples should be thawed immediately prior to the lysis step. Once extracted, RNA should be immediately placed on ice and stored below -20 °C.</p>

Notes

[illegible]

Safety

Reagent	Risk Symbols	Risk and Safety Phrases
Lysis Buffer LB-V		Contains <i>guanidine hydrochloride</i> <ul style="list-style-type: none">• Harmful if swallowed.• Irritating to eyes and skin.• Do not breathe vapors and fumes• Wear suitable protective clothing and gloves.• In case of contact with eyes immediately wash with plenty of water and seek medical advice.• If swallowed, seek medical advice and contact poison control center.

This product is intended solely for research use only. Full material safety data sheet document is available at www.lunanano.com.

For further information about this and other products please visit our website at www.lunanano.com.